

CAVIWIPES for Low Level Cleaning & Disinfecting

****Please use as an addition to the [VCH: Low Level Cleaning & Disinfecting Handout](#)***

CaviWipes:

- Active ingredient: Disobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride
- **Caution:** Corrosive. Avoid contact with eyes, and **MUST** wear gloves when using.
- Wet contact time: 3 minutes
- NOTE: CaviWipes are popular for use in community settings because they are available in individual packs, for convenience for staff. Although this is not the primary VCH approved product, it is still effective against a variety of pathogens.



For more information:

- [VCH Equipment Cleaning Manual](#) for information on products to use on specific items
- [VCH: Low Level Cleaning & Disinfecting Handout](#)