

COVID-19 Doffing PPE for AGMP

Gloves

Remember, the outside of the gloves are contaminated. Grasp palm area of one gloved hand and peel off first glove. Slide fingers of hand under other glove at wrist and peel off. Discard in regular waste¹.

¹If double gloving, remove both pairs.

Hand Hygiene

Clean all surfaces of hands and wrists.

Gown

Unfasten ties (use buddy if available). Pull gown away from the neck and shoulders, touching only the inside of the gown. Turn gown inside out and roll into a bundle. Place in soiled linen hamper. **DO NOT DISCARD IN WASTE UNLESS USING A SINGLE USE GOWN.**

Hand Hygiene

EXIT ROOM

Hand Hygiene

Bouffant Cap

If using, remove by pinching the top and pulling back and away from your face. Do not let the outside surface of the bouffant touch your head while doffing.

Hand Hygiene

Eye Protection

Remove and discard or place in a designated location prior to cleaning and disinfecting.

Hand Hygiene

N95 Respirator

Bend over slightly so that your face is facing towards the floor. Grasp strap from nape of neck and bring it to the front. Grasp strap from crown of head and bring it to the front. Let mask fall away from the face and discard. Do not touch the mask with your hands, only the straps, and do not let the outside of the mask touch your face. Reusable respirators (elastomeric) to be placed in a designated receptacle for reprocessing on removal.

Hand Hygiene

