

COVID-19 Donning PPE for AGMP

Hand Hygiene

Clean all surfaces of hands and wrists.

Gown

Cover torso and wrap around back, fasten in back of neck and waist.

N95 Respirator

Don mask by pulling the lower strap to the nape of the neck and then pulling the upper strap to the crown of the head. Adjust and perform seal check.

Eye Protection

Don either goggles or face shield (the type of face shield that has a strap at the crown of the head, NOT the type attached to a surgical mask).

Bouffant Cap

Optional to control hair outside of OR environments.

Gloves

Extend to cover wrist of gown¹.

¹If double gloving, ensure outer pair covers the wrist of the gown.

