

Infection Prevention and Control

Best Practice Guideline		Patient/Resident Room Cleaning & Disinfection Before UVC Therapy (RD Machine)					
Date		December 9, 2020					
Date Revised							
D	ate Reviewed						
A PRINTED copy of this guideline may not be the most recent version. The OFFICIAL version is located at www.ipac.vch.ca In preparation for patient/resident room cleaning & disinfection, the following must be done:							
1. Clinical Unit staff							
	Remove all patient/resident personal care items (e.g. labelled hospital/LTC supplied liquid soap, mouthwash, tissues, incontinent underwear, etc.). Hospital/LTC supplied personal care items either go home with the patient/resident or are discarded.						
	Discard all leftover supplies that cannot be cleaned and disinfected (e.g. supplies with paper-backed packaging cannot be cleaned and disinfected and must be discarded).						
	Remove mobile medical equipment from room that will not be cleaned & disinfected within the room (e.g. infusion pumps); clean and disinfect 'high touch' points prior to taking to soiled utility room for complete cleaning and disinfecting.						
2.	EVS staff						
	Remove all fabric materials from room (e.g. bedding, pillow cases, towels, gowns, PJs, privacy curtains, etc.). Note: plastic covered pillows (with cloth pillow cases removed) can remain in room.						
	Remove plastic bags from waste containers, linen hampers, and biohazard and chemo waste containers.						
	Systematically clean & disinfect room(i.e. cleanest to dirtiest, top to bottom, left to right or right to left, bathroom is last).						
	Complete 2-step cleaning and disinfection of room:						

December 9, 2020





Infection Prevention and Control

- Step 1: **Clean** equipment/surfaces/rooms with a new wipe(s) or Microfibre cloth(s) using friction (rub/scrub motion) to remove any foreign matter (i.e. dust, soil, food, feces, blood, sputum) and then immediately go to Step 2.
- Step 2: **Disinfect** all equipment/surfaces/rooms using friction (rub/scrub motion) with another new wipe(s) or Microfibre cloth(s).
 - a. To complete the disinfection process, items must be wet enough to maintain the product label stated 'wet contact' or 'dwell time' (e.g. stays wet for 1 minute or 5 minutes). Allow item to <u>air-dry</u> to complete disinfection process. Do not wipe dry to speed things up as this re-contaminates items.
 - b. EVS staff use 4 folds of each of the Microfibre cloths. Each fold is considered to be a clean cloth. Each cloth with 4 folds can clean and disinfect 2 items.

Ц	Expose 'high touch' surfaces as much as possible (e.g. monitors, telephone handles, bed control pad, BP cuff, ceiling lift handles, remotes, etc.) to for allow optimal exposure to UV light.
	Move all furniture away from walls.
	Elevate bed, lift head and foot, make sure bed rails are up and open bed control panel cover.
	Elevate over-bed table.
	Partially open closet and cabinet drawers & doors.

Note: Do <u>NOT</u> place additional items in room (e.g. PPE Cart, Crash Cart) as this will compromise the UVC disinfection.

Items that may be left in room for UVC Disinfection

Bed & Mattress	Bedside Night Stand	Chairs	Over-bed Table	Remote Controls	
Blood Pressure Cuff	Stethoscope	IV Poles	Plastic covered Pi	l Pillows (with linen pillow	
			cases removed)		

