

Vancouver Infection Prevention and Control CoastalHealth VCH Approved Cleaning and Disinfecting Wipes







Product Name	Accel INTERVention™	Accel Rescue™	Ultra Swipes™
Use	Horizontal surfaces, medical devices and equipment	Horizontal surfaces, medical devices and equipment if patient on Contact Plus Precautions ONLY	Medical devices, equipment and electronics where cleaning and disinfection with accelerated hydrogen peroxide is contraindicated
Wet Contact Time	1 Minute 5 Minutes – <i>Candida auris</i>	10 Minutes	1 Minute
PPE	Gloves are recommended for prolonged contact with the product	Gloves	Gloves